

INSTRUCTION MANUAL
FOR
MODEL 4501*
4" FIXED GAP ELECTROMAGNET

CUSTOMER: _____

CUSTOMER ORDER NO: 46700-219621

ALPHA SALES ORDER NO: 060-4180

ALPHA SCIENTIFIC, INC.
460 Roland Way
Oakland, California 94621
Phone: (415) 635-2700

a subsidiary of syston-donner corporation

INDEX

<u>SECTION</u>	<u>TITLE</u>	<u>PAGE NUMBER</u>
I	INTRODUCTION	1
II	GENERAL DESCRIPTION	2 - 3
III	INSTALLATION	4
IV	OPERATION AND PRECAUTIONS	5
V	ELECTRICAL AND COOLING RATINGS	6
VI	MAGNETIZATION CURVES	7
VII	MAINTENANCE	8
VIII	FACTORY SERVICE	9
IX	ACCESSORIES	10
X	4" MAGNET BROCHURE	11

I. INTRODUCTION

The Alpha Model 4501 is a 4" diameter pole electromagnet. The air gap width without pole caps is 3.89 inches. Any desired configuration pole caps are available (See Section IX, ACCESSORIES).

This electromagnet is designed to be used with the Alpha Model 3000 Power Supply (with optional Model 3001 regulator), for experiments such as NMR, EPR and Magnetic Susceptibility measurements, and Hall Effect, Magneto-Optical, Biological studies, etc.

This manual is prepared to acquaint the user with proper installation, operation, maintenance techniques, and to list factory services and available accessories.

II. GENERAL DESCRIPTION

- a. The Model 4501 is intended for use in the laboratory or classroom as a general purpose electromagnet. The magnet, which is rigidly mounted on 45° legs, may be installed on a turntable (Model 4521) for 360° rotation about its vertical axis. When mounted on a turntable, a vernier dial allows rotation to be set within 0.1°. In addition, either magnet or turntable may be installed on a rolling base (Model 4522).
- b. Performance ratings with various pole cap configurations are given on "Magnetization Curves" (Section VI), which show field as a function of air gap and energizing current.
- c. Construction Features
 1. Mechanical - Extreme care has been exercised in the design of this magnet to provide optimum performance. The magnet consists of an "H" type yoke with separate pole bodies and optionally available with water cooled or air cooled coils. The yoke and pole bodies are fabricated from high quality, low carbon steel. The pole caps are made from a highly pure magnetic iron and have been machined and ground to provide maximum field homogeneity. Each pole cap is hard chrome plated to resist wear and corrosion. Using a NMR detector,

the pole caps are shimmed for the most uniform field centered on the pole axis. Axial holes are provided through the pole cores.

2. Electrical - The coils are wound with strip aluminum conductor insulated with mylar and epoxy resin. In the water cooled model, the water cooled heat sink is mounted flat against coil windings for good heat transfer. They are isolated electrically from the windings. The four cooling coils are connected in parallel between common manifolds. The manifolds are terminated with 3/8" female pipe threads.

III. INSTALLATION

- a. If visible signs of damage to crate or to magnet are noticed, call a representative of the transportation company before proceeding further. The magnet must be lifted off the packing skid. Use the lifting eyebolt provided.
- b. The overall dimensions of the Model 4501 are shown in the 4" Magnet Brochure provided. The weight of the magnet is approximately 230 pounds. When the magnet is purchased as a system consisting of both magnet and power supply, two 15' D.C. power supply cables will be provided. When the magnet is purchased alone, D.C. cables are not provided.


IV. OPERATION AND PRECAUTIONS

1. When energizing the magnet, apply full current and then reduce to the operating point as required.
2. To rotate the magnet about vertical axis, unlock the brake before turning the handle.
3. Make sure that the cooling water is turned on.
4. Be sure that there are no loose objects around the magnet while power is on.
5. The magnet cables and hoses should be tightened so that they can not be accidentally broken or disconnected while power is on.
6. Care must be exercised not to damage power cables and interconnecting hoses when rotating the magnet about the vertical axis.

V. ELECTRICAL AND COOLING RATINGS

1. Rated Power = 1 Kilowatt.
2. Rated Voltage = 100 V.D.C. at 10 Amperes.
3. Cooling = 1 G.P.M. at 30 P.S.I.
4. Temperature Rise = 10°C average on water cooled models.
40°C average on air cooled models (at 250 watt)
5. Total Resistance at 20°C = 10 ohms. (Coils in series)


VI. MAGNETIZATION CURVES



HOMOGENEITY PLOT (4.0 K-G.)
4" CYLINDRICAL POLE CAPS ON
MODEL 4501 (WATER COOLED)
FULL SCALE
GAP = 1.50 INCHES

TYPICAL FIELD HOMOGENEITY PLOTS

MODEL 4501 (*)



B = 4000 Gauss at 1" air gap.

MODEL NO. 4501 (WITH MDL 3002-1 P.S.)

4072

FIELD STRAINING: 6.6

20

1 1/2"

RATED POWER: 50


20

TESTS WITH CO.
PARALLEL

SOLE TIP STYLE: CYLINDRICAL

COOLING: WATER (1 G.P.M. @ 30 P.S.I.)

TEMPERATURE RISE: 10°C



MAGNETIZATION CURVE

MODEL NO. 450/1 (*) (WITH MDL 3000 P.S.)

ORDER NO. _____

FIELD STRENGTH: _____ K-GAUSS @ _____ AMPS, WITH _____ GA


RATED POWER: _____ V. D. C. @ _____ AMPS.

POLE TIP STYLE: 4" ϕ CYL.

COOLING: _____

TEMPERATURE RISE: _____


FIELD STRENGTH (KILOGAUSS)



MODEL 4501 (*) -

FIELD STRENGTH - (KILOGAUSS)

AIR GAP - (INCHES)
AT 100 ENERGIZING POWER



VII. MAINTENANCE

- a. Electrical - No special routine maintenance is required.

If the coil insulation should be damaged, check with factory for information regarding repair.

- b. Cooling System - When ordinary tap water is used for cooling, it is recommended that the cooling channels be flushed out periodically (at least annually) with a commercial flushing solution, as used in heat exchangers, to remove slag and deposits.

When the electromagnet is to be stored where freezing is a possibility, flush out the cooling water with compressed air and plug the fittings, or add anti-freeze.

- c. Mechanical

1. To remove pole caps:

- a. Degauss the magnet so that the residual field in the magnet is less than 3 gauss.
- b. Brace the pole cap so that it will not fall when loosened.
- c. Unscrew hex head cap screw holding pole cap to pole.
- d. Carefully remove the pole cap. NOTE: The pole cap may adhere to the pole surface due to the fine surface finish on each part. In this case, it is necessary to slide the pole cap across the pole. The pole cap may be tapped with a soft nylon hammer.

2. To replace pole caps: Reverse above procedure.

VIII. FACTORY SERVICE

1. In all correspondence concerning this electromagnet, please note Model and Serial numbers which are stamped on the specification plate.
2. When requesting replacement parts or ordering accessories, please specify:
 - a. Item or Drawing Numbers, if any.
 - b. Description.
 - c. Number required.
3. Factory service personnel can be contacted by telephone without obligation concerning maintenance or special problems. Maintenance can be performed in the field for a nominal charge. Contact Alpha Sales and Service Department for rates.

IX. ACCESSORIES

Contact the factory Sales Department or your local manufacturer's representative for current price and delivery of accessories.

NOTE: When ordering pole caps, always specify gap desired.

1. Model 4510 4" Diameter Cylindrical Pole Caps
2. Model 4511 4" Diameter Tapered Pole Caps (specify face di
3. Model 4512 4" Diameter Cylindrical Axial Hole Pole Caps
4. Model 4513 4" Diameter Tapered Axial Hole Pole Caps
5. Model 4515 4" Diameter Constant Force $H \, dh/dx$ Pole Caps
6. Model 4516 4" Diameter Constant Gradient dh/dx Pole Caps
7. Model 4520 Variable Height Base
8. Model 4521 Rotating Base
9. Model 4522 Rolling Base
10. Model 1403 Equipment Mounting Table
11. Model 4523 Mobile Cabinet
12. Others Upon Request

Scientific

SUBSIDIARY
BYSTRON  DONNER
CORPORATION


A NEW AND COMPLETE LINE OF 4" MAGNETS

PRECISION LABORATORY
ELECTROMAGNETS
FOR EVERY


Alpha Scientific is now offering a complete line of laboratory electromagnets to meet a wide variety of applications in both industry and research.

A complete family of 4-inch Electromagnets are being offered which are ideally suited to applications such as Magnetic Susceptibility measurements, Electron Paramagnetic, and Nuclear Magnetic resonance demonstrations, and Hall effect studies.


All of the 4-inch electromagnets are designed to be used with Alpha's complete line of Magnetic Susceptibility Systems. Matching Power Supplies and a wide variety of accessories can be purchased for each of these 4-inch electromagnets.



Model 4501 with
Model 4521 Rotating Base
and Model 4522 Rolling Base




Model 4600



Model 4800

MODEL 4501




Features

- * Fixed Air Gap
- * "H" Closed Yoke Design
- * 45° Mounting
- * Water Cooled Coils
- * Easily Removable Pole Caps

Specifications

Power Requirements 1 Kilowatt Maximum
Energizing Current 10 Amperes
Cooling Required 1.0 GPM at 30 psid
Temperature Rise 10°C
Magnet Resistance . . 10.0 Ohms (coils in series)

MODEL 4600



Features

- * Continuously Adjustable Air Gap
- * "H" Closed Yoke Design
- * 45° Mounting
- * Water Cooled Coils
- * Easily Removable Pole Caps
- * Pole Adjustment to Optimize Sensitivity

Specifications


Power Requirements 1.5 Kilowatts Maximum
Energizing Current 30 Amperes
Cooling Required 1.0 GPM at 30 psid
Temperature Rise 15°C
Magnet Resistance . . 1.5 Ohms (coils in series)



MODEL 4800

Features

- * Continuously Adjustable Air Gap
- * "C" or Open Yoke Design
- * Vertical Mounting
- * Water-cooled Coils
- * Removable Pole Shoes



Specifications

Power Requirements.....2 Kilowatts Maximum
 Energizing Current.....50 Amperes
 Cooling Required.....1.0 GPM at 30 psid
 Temperature Rise.....20°C at Full Power
 Magnet Resistance...0.8 Ohms (coils in series)

Physical Characteristics






Air Gap.....0 to 4.25
 Pole Diameter.....4.0" Nominal
 Coil Spacing.....6⁵/₈"
 Axial Hole.....2⁵/₃₂" Diameter (thru Pole Body)
 Dimensions.....Height, 19⁵/₈"; Width, 28³/₁₆";
 Depth, 16⁷/₈"

Centerline Height:


- 11¹/₈" with Vertical Mounting
- 17¹/₈" Mounted on Rotating Base
- 19⁷/₈" Mounted on Rotating and Rolling Base

Shipping Weight:

- 585 lbs. with Vertical Mounting
- 685 lbs. with Rotating Base
- 715 lbs. with Rotating and Rolling Base



B=3500 Gauss at 3/4" air gap.



B=6800 Gauss at 3/4" air gap.

TYPICAL FIELD HOMOGENEITY PLOTS
 TAKEN ON MIDPLANE
 SCALE REDUCED TO 1/4

Matching Power Supplies For High Laboratory Magnets

Magnet	Power Supply	
Model 4501.....	Model 3000, 100 volt, 10 Amp Filtered, Unregulated	Shipping Weight: 80 lbs.
	Model 3001, Regulator (5×10^{-4})	Shipping Weight: 30 lbs.
Model 4600.....	Model 3002-1, 50 volt, 30 Amp Regulated to $\pm 5 \times 10^{-5}$	Shipping Weight: 170 lbs.
Model 4800.....	Model 3003-1, 40 volt, 50 Amp Regulated to $\pm 1 \times 10^{-5}$	Shipping Weight: 250 lbs.

Accessories

- Model 4510 Cylindrical Pole Caps
- Model 4511 Tapered Pole Caps
- Model 4512 Cylindrical Axial Hole Pole Caps
- Model 4513 Tapered Axial Hole Pole Caps
- Model 4514 Pole Shim Pairs
- Model 4515 Constant H dH/dx Pole Caps
- Model 4516 Constant dH/dx Pole Caps
- Model 4520 Variable Height Base
- Model 4521 Rotating Base
- Model 4522 Rolling Base
- Model 1403 Equipment Mounting Table
- Model 4523 Mobile Cabinet

Ordering Information

- Please Specify:
- Model Numbers of units desired
 - Magnet Price includes one pair of cylindrical or tapered pole caps (specify)
 - Pole Cap requirements (Face Dia. Gaps., etc.)
 - Correct shipping address
 - Method of shipment
 - Special instruction if any



C O R P O R A T I O N

460 ROLAND WAY, OAKLAND, CALIFORNIA 94621 (415) 635-2700